

Lancing Device & Lancets





Universal

Universal Design for fit with most common lancets

Ultra-thin, electro-polished lancets for comfortable sampling

EASY

- Low fore to load the device with convenient end grip
- Adjustable depth setting easy to rotate and with clear visual marking

GENTLE

- Compact size with tight linear lancet guiding for gentle puncture
- Quite and soft activation for seamless blood sampling

SAFE

- Color-Coded loading indicator to distinguish when the device is ready to be used
- Automated ejector to safely remove lancet from the device

UNIVERSAL & ULTRA-FINO

- 23G to 33G ultra-thin lancets of universal fit
- Tri-Bevel, electro-polished lancets for enhanced user comfort



Automated Ejector for safe lancet removal



INSTRUCTIONS FOR USE

■ Intuitive and easy to follow steps for user comfort.



Remove the Device Cap



Insert a new Lancet Into the Device



Twist off and pull out safety cap of the Lancet



Put the Device Cap back on



Adjust the lancing depth



Load the Device (Click Indicator)



Activate the Device to obtain blood sample



Drive the Lancet into the safety cap



Eject the lancet safetly from the device

ORDERING INFORMATION

CODE	DIAMETER (G) / DEPTH (mm/in)	DESCRIPTION	PACKAGING
1200	10 - Depth Setting	Universal Lancing Device	1 pc/ box 100 pcs / case
1201	23G	23G Universal Blood Lancets	100 pcs / box 10.000 pcs / case
1202	26G	26G Universal Blood Lancets	100 pcs / box 10.000 pcs / case
1203	28G	28G Universal Blood Lancets	100 pcs / box 10.000 pcs / case
1204	30G	30G Universal Blood Lancets	100 pcs / box 10.000 pcs / case
1205	31G	23G Universal Blood Lancets	100 pcs / box 10.000 pcs / case
1206	32G	32G Universal Blood Lancets	100 pcs / box 10.000 pcs / case
1207	33G	33G Universal Blood Lancets	100 pcs / box 10.000 pcs / case

WWW.ONE-CARE.COM

