

ONE-CARE[®]

SAFETY PEN NEEDLES **PRO DUO**



Your Ideal Choice for Duo Safety Pen Needles!



Safe

Improved visibility during priming, automatic shielding and red color indicator after use

Universal


Compatible fit with most pen injectors

Dual Safety

Back-end needle automatically covered after removal from pen

Ergonomic

Ergonomic main housing for comfortable grip



SAFE

- Needle automatically shielded before and after use to prevent accidental needle-sticks
- Low resistance during injection and activation of the safety mechanism
- Red-color activation indicator to differentiate before used and non-used device


CONVENIENT

- See-through window for improved visibility during priming
- Ergonomic design for easy attachment into pens
- Optimal skin contact platform that promotes 90-degree injection angle

UNIVERSAL

- Compatible with most popular pens on the market
- Different needle configurations available


WWW.ONE-CARE.COM



ONE-CARE[®]
SAFETY PEN NEEDLES **PRO DUO**

INSTRUCTIONS FOR USE

- Easy-to-use design with intuitive handling and simple steps





COMPATIBILITY

Sanofi Aventis	ClikSTAR, SoloSTAR, JuniorSTAR			
Lilly	KwikPen, HumaPen Memoir, HumaPen Luxura HD, HumaPen Luxura, HumaPen Savvio, HumaPen Ergo II			
Novo Nordisk	FlexTouch, NovoPen 3, 4 and 5			
	Innolet, FlexPen			
	NovoPen Echo, Victoza			
Novo Nordisk	Autopen Classic, Autopen 24			
Berlin Chemie	BerliPen Areo, BerliPen Junior			

ORDERING INFORMATION

- Available in different needle configurations :

REF	NEEDLE	COLOR	PACKAGING
3060	ONE-CARE PRO DUO Safety Pen Needles 31G x 5mm		100 pcs / box 1200 pcs / case
3061	ONE-CARE PRO DUO Safety Pen Needles 31G x 8mm		100 pcs / box 1200 pcs / case

WWW.ONE-CARE.COM

MediVena[®]
inspired by innovation

www.medivena.com

AMERICAS:
MediVena
1 Glenlake Parkway, Ste 700, Atlanta, GA 30328, USA
+1 404 514 2586 info@medivena.com
CUSTOMER SERVICE TOLL FREE: 888 270 6457

EUROPE:
MediVena Sp. z o.o.
Al. Jana Pawła II 27, 00-867 Warsaw, POLAND
+48 695 995 700

